



WATERCARE

Reverse Osmosis Restart Guide

Reverse Osmosis Units (RO)

ROBustPro, ROBust1000, ROBust1500, ROBust3000:

Flush water through the unit for between 30 – 60 minutes or until 50 litres have been flushed.

To do this please follow the below steps:

1. Ensure the isolation valve on the outlet line is in the closed position.
2. Remove the push fit line from the outlet of the isolation valve.
3. Position the outlet isolation valve into a sink, waste pipe, or bucket.
4. Open the isolation valve and allow it to run for approximately 30 – 60 minutes or to have reached 50 litres in output water (whichever is soonest).
5. Close the isolation valve, and reconnect to the outlet supply line.



If you need any help flushing your RO unit, please call us on **01279 780250** (opt. 2) or email: service@watercare.co.uk.