



WATERCARE

WHAT TO DO IF YOU HAVE NOT USED YOUR REVERSE OSMOSIS WATER TREATMENT UNIT FOR OVER 4 WEEKS

Important hygiene advice

Reverse Osmosis Units (RO)

ROBustPro, ROBust1000, ROBust1500, ROBust3000:

Flush water through the unit for between 30 – 60 minutes or until 50 litres have been flushed.

To do this please follow the below steps:

- Ensure the isolation valve on the outlet line is in the closed position.
- Remove the push fit line from the outlet of the isolation valve.
- Position the outlet isolation valve into a sink, waste pipe, or bucket.
- Open the isolation valve and allow it to run for approximately 30 – 60 minutes or to have reached 50 litres in output water (whichever is soonest).
- Close the isolation valve, and reconnect to the outlet supply line.



If you need any help flushing your RO Unit, please call us on
01279 780250
(opt. 2) or email:
service@watercare.co.uk